

aureliagloves.com
1.877.AURELIA

AURELIA® GLOVES

  
@aureliagloves



WHEN YOU NEED HEAVY DUTY
PERFORMANCE, THERE'S ONLY ONE
GLOVE THAT GETS THE JOB DONE.

IGNITE THE WORKDAY



TEXTURED



7.0 mil



AMBIDEXTROUS



HEAVY-DUTY



LATEX FREE

IGNITE™

NITRILE GLOVES

- 7 Mil Orange Powder Free Nitrile Examination Gloves
- Ambidextrous with “Max Grip” Texture
- Exceptional Tactile Sensitivity and Comfort
- For Single Use Only
- Great for EMT, Industrial and Automotive Applications



DIMENSIONS

SIZE	GLOVE WEIGHT (g)	PALM WIDTH (mm)		LENGTH (mm)	
		ASTM D6319	IGNITE®	ASTM D6319	IGNITE®
M	8.4±0.2	95±10	90-99	Min. 230	Min. 240
L	9.0±0.2	110±10	100-109	Min. 230	Min. 240
XL	9.8±0.2	120±10	110-119	Min. 230	Min. 240
XXL	10.5±0.2	-	120-129	-	Min. 240

THICKNESS

AREA OF THICKNESS MEASUREMENT	SINGLE WALL (mils)	
	ASTM D6319	IGNITE®
Finger	Min 0.05	Min 0.10
Palm	Min 0.05	Min 0.8
Cuff	N/A	Min 3.6

PHYSICAL PROPERTIES

	BEFORE AGING		AFTER AGING	
	ASTM D6319	IGNITE®	ASTM D6319	IGNITE®
FINGER	Min 14	Min 14	Min 14	Min 14
PALM	Min 500	Min 500	Min 400	Min 400
CUFF	N/A	Min 9	N/A	Min 6

AURELIA® GLOVES

aureliagloves.com | 1.877.AURELIA



@aureliagloves