

Bold

Aurelia® Bold® Black Nitrile
Powder Free Examination Gloves



**Durable Black Nitrile
Maximum Strength**

- 5 mil thickness for improved puncture resistance
- Finger-textured finish for improved gripping action
- 100% latex free, black color




AURELIA®
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Bold

Product Description

Aurelia® Bold® Black Nitrile (Non Latex) Powder Free Examination Gloves are made from 100% Nitrile (Non Latex), ambidextrous, black color, Finger Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Bold® Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	26.1	550
Small	26.1	550
Medium	26.5	550
Large	26.1	595
Extra Large	26.2	580
ASTM D6319-00	14 MPa minimum	500% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	24.0	550
Small	24.0	550
Medium	24.0	555
Large	24.0	585
Extra Large	23.7	570
ASTM D6319-00	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	4.5 / 0.11 3.2 / 0.08 minimum	5 / 0.12 3.2 / 0.08 minimum
Small Aurelia Bold ASTM D6319-00	9.4 / 240 8.7 / 220 minimum	4.5 / 0.11 3.2 / 0.08 minimum	5 / 0.12 3.2 / 0.08 minimum
Medium Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.11 3.2 / 0.08 minimum	5 / 0.12 3.2 / 0.08 minimum
Large Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.11 3.2 / 0.08 minimum	5 / 0.12 3.2 / 0.08 minimum
Extra Large Aurelia Bold ASTM D6319-00	9.4 / 240 9.1 / 230 minimum	4.5 / 0.11 3.2 / 0.08 minimum	5 / 0.12 3.2 / 0.08 minimum

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

Extra Small	Small	Medium	Large	Extra Large
73995	73996	73997	73998	73999

Authorized Aurelia Distributor:



AURELIA Div. Supermax Inc.

1899 Sequoia Drive, Aurora, IL 60506 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
Web: www.aureliagloves.com | Email: info@aureliagloves.com | Facebook: www.facebook.com/aureliagloves