

# Distinct

Aurelia® Distinct™ Honeycomb® Textured  
Powder Free Latex Examination Gloves



**Unique Honey Comb Texture  
Superior Tactile**



- Patented Honeycomb texture for extraordinary gripping
- Polymer-coating enhances donning and keeps hands smooth
- Excellent tactile sensitivity



**AURELIA®**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Distinct

## Product Description

Aurelia® Distinct™ Honeycomb® Textured Powder Free Latex Examination Gloves are made from 100% natural Malaysian rubber, ambidextrous, creamy white color, Honeycomb® Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Distinct™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

## DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Small</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Medium</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Large</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Extra Large</b> Aurelia Distinct ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

### ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
29225	29226	29227	29228	29229

Authorized Aurelia Distributor:

## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



**CAUTION STATEMENT:** This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

**AURELIA Div. Supermax Inc.**

2225 White Oak Circle, Aurora, IL 60502 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
 Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aueliagloves.com](mailto:info@aueliagloves.com)