

Refresh


Aurelia® Refresh™ Honeycomb® Textured
Peppermint Powder Free Latex Examination Gloves



Mint Scented
Exclusive Honeycomb Texture



- Refreshing peppermint scent for you and your patient
- Patented Honeycomb texture for extraordinary gripping and tactile feel
- Polymer-coated enhances donning and keeps hands smooth



AURELIA
Sign of Comfort™

www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Refresh

Product Description

Aurelia® Refresh™ Honeycomb® Textured Powder Free Latex Examination Gloves are made from 100% Natural Malaysian rubber, ambidextrous, green color with peppermint scent, Honeycomb® Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Refresh™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
Extra Small Aurelia Refresh ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Small Aurelia Refresh ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Medium Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Large Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
Extra Large Aurelia Refresh ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
99225	99226	99227	99228	99229

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor:

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

AURELIA Div. Supermax Inc.

2225 White Oak Circle, Aurora, IL 60502 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
Web: www.aureliagloves.com | Email: info@aureliagloves.com